

Certification of compliance

Product : Milk Twinsensor BT Demo Kit
Simultaneous Detection of Betalactam & Tetracycline residues in Milk
 Référence : BT00640
 Number of tests : 48 (Demo Kit)
 Method : 3+3 min 50°C

Batch : 050305
 Exp date : 05NOV2005

approved QC

First instantaneous reading

Ampicillin ng/ml or PPB	Visual interpretation
0	negative
1	negative
2	negative
3	Low positive
4	positive
5	positive
6	positive
8	positive
10	positive
20	positive

Tetracycline ng/ml or PPB	Visual interpretation
0	negative
10	negative
20	negative
30	negative
40	Low positive
50	positive
60	positive
80	positive
100	positive
200	positive

Antibiotics used : WHO/APL : <http://www.apl.apoteket.se/Engelska/who/index.htm>
 - Ampicillin sodium : cat. Number 9930212
 - Tetracycline Hydrochloride : cat. Number 9930429
 - Fresh Raw Milk.

Note : Nothing to report

QUALITY ASSURANCE

Benoît LEMMENS

Date

5MAY2005 unisensor s.a.
 Liège Science Park
 B-4031 Angleur BELGIUM